

Practical:

Single crystals and powders

Agenda

- ★ *Laue Camera*
- ★ *Use and modify the diffractometer*
 - *PSI_DMC*

2019 CSNS

McStas
School

McStas



Laue Camera

First insert a source

File → New Instrument

Insert → Sources → Source_simple

// insert components here (e.g. Insert -> Source -> ...)

```
COMPONENT source_simple = Source_simple(
    radius=0.05, dist=5, focus_xw=0.02,
    focus_yh=0.05, lambda0=2, dlambda=1.9)
```

```
AT (0, 0, 0) RELATIVE PREVIOUS
```

2019 CSNS

McStas

School

McStas



2019 CSNS

McStas

School

McStas



Laue Camera

Now add a guide

Insert → Optics → Guide

```
COMPONENT guide = Guide(
    w1=0.02, h1=0.05, w2=0.02, h2=0.05, l=20, m=1)
AT(0,0,2) RELATIVE source_simple
```



Laue Camera

Add a sample – in this case a standard crystal

Insert → Samples → Single_crystal



```
COMPONENT single_crystal = Single_crystal(
    reflections="Al2O3_sapphire.lau",
    yheight=0.05, radius=0.01, mosaic=1, delta_d_d=1e-4,
    az=4.757, ay=0, az=0, bx=2.3785, by=0, bz=-3.364,
    cx=0, cy=12.9877, cz=0,
    p_transmit=0.1, order=1)
```

AT (0, 0, 20.5) RELATIVE PREVIOUS



Laue Camera

Add *the* ideal Laue Camera Monitor – covering 4π



Insert → Monitors → PSD_monitor_4PI

```
COMPONENT fourpi = PSD_monitor_4PI(  
    radius=1, filename="fourpi.dat", nx=201, ny=201)  
AT (0,0,0) RELATIVE PREVIOUS
```

2019 CSNS


McStas
School

McStas




Laue Camera

Run your simulation (you can safely increase the number of rays to 10^7 - the ncount) – you should get something like:



Laue Camera



2019 CSNS
McStas
School

McStas



The coherent scattering is not much stronger than the incoherent “background”. Let’s use EXTEND and WHEN to make a monitor which only senses the coherent signal.

- In the DECLARE section of your instrument – declare a variable of **type char**;
- Add an EXTEND-block to the end of the Single_crystal component:
- **EXTEND**
- `% {`
- `myvar = hkl_info.type;`
- `% }`
- Add another 4-PI monitor and insert **WHEN** (myvar==99) in front of the AT keyword. (99 is the ascii-code for ‘c’)
- Think of a possible way to only monitor incoherent scattering... HINT: type=='i' for incoherent

Laue Camera

- ▶ Play around with this example instrument:
 - ▶ Add an `Arm` components before the sample to allow rotation around the Y-axis.
 - ▶ Add "`SPLIT 20`" before the sample `COMPONENT` statement. What happens?
 - ▶ Add a set of arms before the sample to add Y, Z, Y rotations (Eulerian cradle). Make the rotation angles input parameters.
 - ▶ Insert a different crystal instead e.g. "`Al.lau`". i.e. change the crystal unit cell parameters and the reflection list.
 - ▶ Insert a powder sample instead

PowderN and More



中国散裂中子源

2019 CSNS
McStas
School

McStas



- In this exercise we will try to put two powder samples together in a few ways and compare the results.*
- 1) Two samples in the same spot, and stochastically choose between them*
- 2) Two samples stacked vertically*
- 3) Two samples where one is behind the other.*


中国散裂中子源

2019 CSNS
McStas
School

McStas



PowdersN and more



SCHOOL

McStas



- Let's use the *PSI_DMC* instrument as a starting point. We will now make the simulation randomly choose between two powders. The instrument file can be found in the McStas distribution. (File → New from template → PSI → PSI_DMC)

- 1) Add another powder in the same spot as the one already there. Look in the mcstas data directory for sample .laz-files
- 2) Add **double r;** inside the **DECLARE** section of the instrument file.
- 3) Add an Arm in front of the first one, and add to it an **EXTEND**-block. Add the following code in it: **r=rand01 () ;**
- 4) Now add the following before the **AT** on the two powders. **WHEN (r<0.5)** and **WHEN (r>0.5)** respectively.
- 5) Run the instrument – Do you get what you expect?
- 6) What would you change to make the mixing factor !=0.5?
- 7) Use a similar technique to the Laue camera to make the detector only catch scattering from one sample.



PowderN and More

- *Change this to have two samples on top of each other.*
 - 1) *Make a new copy of the instrument (or remove the edits you did before, but leave the second PowderN sample in place).*
 - 2) *Change the y-position and size of the samples to be:
 $\pm \text{sample_height}/4.0$ and $\text{sample_height}/2.0$ respectively*
 - 3) *Add the statement **GROUP samples** after the AT at both samples. (N.b. “samples” is a name chosen arbitrarily. It has to be different than the component names though.)*
 - 4) *Run a simulation – is there any difference to the previous result? Why/Why not?*

PowderN and More

- *Move the samples around such that one is in front of the other.*
 - 1) *Run a simulation – Do you still see the signatures of both samples?*
Do you remember why this can be?
 - 2) *How can we get around this?*



中国散裂中子源

2019 CSNS
McStas
School


McStas



中国散裂中子源

2019 CSNS
McStas
School

McStas




Intermission:

A quick trick to remove the direct beam

- *If your monitor also can be hit by the direct beam, “swamping” the signal, you can do this:*
- *Add the following code just after your sample code:*

```
EXTEND
%{
    if (!SCATTERED) {ABSORB;}
%}
```
- *This will terminate all rays which the sample-code has not flagged as scattered. Note that the McStas definition of scattered includes many things (guide-wall reflections etc.)*

PowderN and More



Increase the height of the detector and make it resolve the signal along y.

Set:

Options="banana, theta y auto limits bins=20", yheight=0.3


...In the detector.

Change the sample to be e.g. an Aluminium crystal.

```
COMPONENT single_crystal = Single_crystal(
    reflections="Al.lau",
    yheight=0.05, radius=0.01, mosaic=1, delta_d_d=1e-4,
    az=4.0495, ay=0, ax=0, bx=4.0495, by=0, bz=0,
    cx=0, cy=4.0495, cz=0,
    p_transmit=0.1)
```

AT (0, 0, 0) RELATIVE PREVIOUS

PowderN and More



Insert more detector banks to modify DMC to become more like a Laue camera (catch more of the crystal signal). New banks could be above and below and on the other side of the sample. You will need to use a **GROUP** for this.

Increase the mosaicity of the crystal to scatter more of the beam. This can be done to an extent – think about what limits it, and how you could extract such limits from the simulations.

Using the WHEN keyword we can make a “sample changer”. Can you think of how?

HINT: similar to the earlier exercise with a mix of two crystals.