




# From modelling to McStas and more

## A short guide

E. Farhi, *ILL*



# Data formats





# McStas Data formats

McStas uses a number of input data formats, all text based.

- **Lau** (Single\_crystal, PowderN, Isotropic\_Sqw)
- **Laz** (PowderN, Isotropic\_Sqw)
- **Sqw** (Isotropic\_Sqw)
- **qSq** (Isotropic\_Sqw) – mostly for liquids
- **nxs** (Sample\_nxs) – this is ***not*** a NeXus/HDF
- **PDB** (SANSPDB)

The question is how to generate these files...




*By hand ?*





# Material modelling formats



- MD (classical or DFT)
  - trajectories, many data formats
- Lattice dynamics:
  - PhonoPy FORCE\_SETS
  - ...



# Nuclear data base formats



- ENDF, JEFF, JENDL, ...
  - all use a standard format for a set of materials. Most are monoatomic / incoherent
- ACE
  - Same as ENDF, but in binary for MCNP
- Ncrystal (2.5)
  - NCMAT, NXS, LAU





# Experiments

- Intensity( $d$ )  $\rightarrow$   $F(Q)$   $\rightarrow$   $S(Q)$  and Laz
- Intensity(HKL)  $\rightarrow$   $F_2(\text{HKL})$   $\rightarrow$  Lau
- Intensity(theta,ToF)  $\rightarrow$   $S(q,w)$
- Often saved as text and/or NeXus








# Tools

- iFit
- Mantid
- Python/ASE (IO)
- MDANSE/nMoldyn (converters)
- Scripts of all kinds
  - OpenBABEL, ...
- cif2hkl
- bash





# Material DataBases

- ICSD <<http://icsd.ill.fr/icsd/index.html>>
- COD <<http://crystallography.net>>
- <<https://www.materialsproject.org/>>
- PhononDB <<http://phonondb.mtl.kyoto-u.ac.jp/>>
- Nomad <<https://repository.nomad-coe.eu/>>

Mostly provide structures...






# McStas: structure: powder



- For powders, the pure elastic diffraction is easy to prepare: LAZ
  - PowderN
  - Single\_crystal(powder=1, ...)
  - Isotropic\_Sqw



CIF file  
FullProf/ShelX

Terminal:  
`cif2hkl --powder file.cif`

McStas Laz file

CIF file  
FullProf/ShelX

iFit:  
`cif2hkl('file.cif', '', 'p')`

McStas Laz file





# McStas: structure: SX

- For single crystals, the pure elastic diffraction is easy to prepare: LAU
  - PowderN
  - Single\_crystal
  - Isotropic\_Sqw

CIF file  
FullProf/ShelX

Terminal:  
`cif2hkl --xtal file.cif`

McStas Lau file

CIF file  
FullProf/ShelX

iFit:  
`cif2hkl('file.cif', '', 'x')`


McStas Lau file





# McStas: structure: COD

- cif2hkl from iFit can input chemical formula
- Requires internet connection to COD
- Hill notation: C H then alpha order



Chemical  
formula


iFit:  
cif2hkl('cod: La Mn O3', '', 'p')

McStas Laz file

Chemical  
formula

iFit:  
cif2hkl('cod: O3 Sr Ti', '', 'x')

McStas Lau file




# McStas: structure: others



- Other McStas structure file formats
  - SANSPDB and similars: get from Data Bases
  - Sample\_nxs: manual edit only (nxs, text)
  - Isotropic\_Sqw: manual edit only (qSq, text)
  - NCrystal\_sample\*: (ncmat,nxs,lau)

\*in McStas from v. 2.5



OpenBabel supports many data files





# McStas: dyn from ToF spec

- Import and correct data with e.g. LAMP, Mantid
- Integrate radially ( $|q|$ ) and compute  $S(q, \omega)$
- Export to e.g. HDF/NeXus or other
- Import then export with iFit:
  - `sqw=iData_Sqw2D('file')`
  - `sqw.saveas('file.sqw', 'mcstas')`
- Can/should always edit manually (text)

```
# weight      4.002    in [g/mol]
# density      0.4784   in [g/cm^3]
# sigma_abs    0.00747  absorption scattering cross section in [barn]
# sigma_coh    1.34     coherent scattering cross section in [barn]
# sigma_inc    0         incoherent scattering cross section in [barn]
# Temperature  2         in [K]
# classical    0         experimental, contains Bose factor
```



# McStas: dyn from exp. DOS

- Import and correct data with e.g. LAMP, Mantid
- Estimate  $g(\omega)$  with e.g. LAMP, Mantid, MUPHOCOR
- Export to e.g. HDF/NeXus or other
- Import then export with iFit:
  - `g=iData_vDOS('file')`
  - `inc=g.incoherent;`
  - `inc.saveas('inc.sqw','mcstas')`
- Can/should always edit manually (text)







# McStas: dyn from nuclear DB

- Import and export with iFit:
  - `sqw=iData_Sqw2D('tsl.endf')`
  - `saveas(sqw(1), 'inc.sqw', 'mcstas')`
- Ncrystal (2.5) will provide an alternate solution
- Can/should always edit manually (text)

```
# weight      4.002    in [g/mol]
# density     0.4784   in [g/cm^3]
# sigma_abs   0.00747  absorption scattering cross section in [barn]
# sigma_coh   1.34     coherent scattering cross section in [barn]
# sigma_inc   0        incoherent scattering cross section in [barn]
# Temperature 2        in [K]
# classical   0        experimental, contains Bose factor
```





# McStas: dyn from MD

- Import trajectory into MDANSE (converters)
- Double click the trajectory to allow analysis
- Use Plugins:
  - Analysis/Scattering/DCSF and DISF
  - *Analysis/Scattering/Structure Factor*
  - *Analysis/Dynamics/Density of States*
- Import then export with iFit:
  - `sqw=iData_Sqw2D('DCSF.nc')`
  - `sqw.saveas('coh.sqw','mcstas')`
- Can/should always edit manually (text)





# McStas: dyn from Lattice Dynamics

- Using PhonoPy:
  - iFit: `s4d=sqw_phonons('phonopy_directory')`
- Computing on the fly:
  - iFit: `s4d=sqw_phonons('cod: La Mn O3')`
  - `p=powder(s4d)`
  - `d=iData_Sqw2D(iData(p, [], 0:10, 0:100))`
  - `d.saveas('coh.sqw','mcstas')`
- Other:
  - Import in iFit, convert to Sqw2D, export.






# Neutron events

- McStas handle ‘rays’ as packets of neutron events.
- Summing-up builds intensity and distributions.
- McStas can read/write ‘events’ from other software:
  - MCNP TRAC and SSW
  - TRIPOLI
  - MCPL\*
- McStas can generate events as well for Mantid via NeXus/HDF files.

\*Allows interchange with many other MC codes: MCNP, Geant4, SIMRES, Vitess, ...




# McStas automation: bash



- Can use cif2hkl and mcrun commands:
  - **cif2hkl** 'file' --out 'file.laz' --powder
  - **mcrun** --dir=test **instr** Powder=file.laz
  - **mcplot.pl** -psc test

```
[~]$ bash
```





# McStas automation: iFit

- iFit is a generic infrastructure which gathers data sets handling, fitting,  $S(q,w)$  models, and McStas hooks.

*iFit*







# McStas iFit: build

- McStas can be controlled from within iFit.
- Open Matlab/iFit
- Create the McStas model with:
  - `model = mccode('instr')`
- Plot the geometry with:
  - `plot(model)` *% has contextual menus*
- Edit the instrument and re-compile
  - `edit(model)`





Figure 2: TextEdit: /tmp/tp109e7fc7\_fd80\_

```
File Edit Help
McStas instrument definition URL=http://www
Instrument: Template monochromator Diffractometer
%Identification
Written by: E. Farhi
Date: 13 Apr 2006
Origin: LLB/ILL
Release: McStas CVS_080624
Version: $Revision$
%INSTRUMENT_SITE: Templates
Simple monochromator Diffractometer for powders
%Description
```

# McStas iFit: eval


- Run with (default 1e6 event and pars)
  - `data = iData(model, [], nan);`
  - `subplot(model) % plot monitors, contextual menus`
- Specify parameters
  - `data = iData(model, 'lambda=2.36; coh=Cu.laz')`
- Do a scan:
  - `data = iData(model, 'lambda=[1.2 2.4 3.6]')`
- Change neutron events #
  - `model.UserData.options.ncount = 1e7;`





# McStas iFit: optim


- Fix all parameters but *lambda*, Maximize model value:
  - `mlock(model, 'all');` `munlock(model, 'lambda')`
  - `xlim(model, 'lambda', [1 3]);` % bounds
  - `fmax(model, 'lambda=2.36', '', nan)`
  - `fmax(model, 'lambda=2.36', 'OutputFcn=fminplot', nan)`





# McStas iFit: advanced

- You can add McStas models:
  - `model=mccode('instr1')+mccode('instr2') + ...`
- For instance:
  - Instr1: structure
  - Instr2: spin-wave
  - Instr3: phonons
  - Instr4: incoherent
- The different models can also be assembled as a set of samples chosen at execution in a single McStas instrument.





# Links to software

---

- McStas <http://www.mcstas.org> (includes cif2hkl)
- iFit <http://ifit.mccode.org> (includes cif2hkl)
- Mantid <http://www.mantidproject.org>
- MDANSE <http://www.mdanse.org>
- ASE <https://wiki.fysik.dtu.dk/ase/>
- Ncrystal <https://github.com/mctools/ncrystal>

