

Virtual Experiments

LaMnO₃ sample: adding all contributions

E. Farhi, *ILL*

Our goal

- Model a ToF with a LaMnO₃ sample
- This material has contributions:
 - Structure (Bragg peaks)
 - Phonons
 - Spin-waves
 - Incoherent scattering
- We shall prepare all this and actually assemble a full model.
- Dropbox: 19th: VE LaMnO₃



Limitations

- We shall use a simple ToF direct spectrometer.
- The sample will be assumed in powder state.
- The intensity of the contributions will have to be scaled with the proper cross sections. Will not be absolute.
- No phonon-spinon coupling as contributions are computed separately.
- Antiferromagnet, with localized coupling (BZ boundary).




iFit: data sets



- Import a data file:
 - `data = iData('filename')` % many formats
 - `data = iData(x,y,z,..., signal)`
- Convert to a flavour:
 - `data = iData_Sqw2D(data)`
- Can apply many operators between data sets
 - `+ - / *` trigo ... 193 operators.



iFit: models



- Create a model:
 - model = iFunc('p(1)*x+p(2)')
 - model = gauss
- Convert to a flavour:
 - model = iFunc_Sqw2D(model)
- Can apply many operators between models
 - +-/ * trigo ... 81 operators.



iFit: predefined models

- Predefined models: 57
 - `[o,m] = fits(iFunc)` % optimizers, models
 - `doc(iFunc, 'Models')`
- Those we shall use here:
 - **sqw_phonons**: $S(q,w)$ 4D lattice dynamics / DFT
 - **sqw_spinwave**: SPINWAVE model




sqw_phonons

- Uses transparently ASE+PhonoPy to compute lattice dynamics estimate of the Forces in a SX.
- Produces a 4D $S(q,w)$ model [iFunc_Sqw4D]
- Supports all ASE calculators, e.g. QE VASP ...
- For complex crystals, computation can be days.
- iFunc Sqw4D flavour methods:
 - vDOS, thermochemistry, plotting, powder average, ...



sqw_phonons (2)

- Will use MPI and all CPU's
- Can input CIF/PDB/... file, as well as common ASE Atom builder (python)
- Can also input a 'cod: formula'
- **Syntax:** `sqw_phonons('crystal', 'calculator', ...)`
 - `sqw_phonons('Al.cif', 'qe')`
 - `sqw_phonons('cod: Al', 'vasp')`
 - `sqw_phonons('bulk("Al", "fcc", a=4.05)', 'gpaw')`



Exercise 1: *sqw_phonons Al*

- Generate the ‘default’ Al $S(q,w)$ model:
 - `Al = sqw_phonons('bulk("Al", "fcc", a=4.05)', 'emt')`
- Plot the dispersion and the BZ cube
 - `plot(Al)` *% dispersion and DOS*
 - `plot3(Al)` *% 3D cube*
 - `slice(Al)` *% play with isosurface and cuts*
- The views have contextual menus on figure, axis and plots.



sqw_spinwave

- Uses transparently SPINWAVE.
- Has a special ‘powder’ mode (that we shall use).
- Can define variable model parameters *\$par*
- Produces a 2D/4D $S(q,w)$ model [iFunc_Sqw2/4D].
 - *sqw_spinwave*(‘file’)
 - *sqw_spinwave*(‘file’, ‘powder’)
 - *sqw_spinwave*(‘file’, ‘edit’) % then insert *\$par*, e.g. J1...




LaMnO₃: Structure

- Just need to generate a LAZY/PulverX Laz file
- Matlab/iFit:
 - `cif2hkl('cod: La Mn O3', 'LaMnO3.laz')`
- Then use it with e.g. **PowderN**




LaMnO₃: Phonons

- Calculation done with QE:
 - `s=sqw_phonons('cod: La Mn O3','qe')` % 15 h, 24 cpus
- Start Matlab/iFit:
 - `load LaMnO3`
 - `disp(s)` % a 4D model
 - `figure; plot3(s)` % use contextual menu 'invert transparency'
 - `slice(s)`



LaMnO₃ → Sqw phonons coherent

- Matlab/iFit:
 - p=**powder**(s) % *powder average model*
 - q=0.01:0.025:4; w=0:0.05:50;
 - coh=**iData**(p,[], q , w);
 - coh=**iData_Sqw2D**(coh);
 - figure; **plot**(coh)
 - **saveas**(coh, 'LaMnO₃_pow_coh.sqw', 'mcstas')




LaMnO₃ → Sqw incoherent

- Matlab/iFit: with same (q,ω) range as coherent.
 - g=`dos`(s, w); % vDOS with coh 'ω' values
 - `plot`(g)
 - inc = g.`incoherent`('q', q, 'm', 242, 'T', 10);
 - figure; `subplot`(inc)
 - sinc = `deBosify`(plus(inc));
 - ylim(sinc,[-0.15 50]); sinc{1}=w; sinc{2}=q;
 - `saveas`(sinc, 'LaMnO₃_pow_inc.sqw', 'mcstas')



LaMnO₃: SpinWaves → Sqw

- Create the spin-wave model:
 - SW=sqw_spinwave('LaMnO₃.txt','powder')
- Evaluate on our axes:
 - dSW=iData(SW, [], q, w); % **SLOW** (20 min)
 - dSW=iData_Sqw2D(dSW);
 - dSW.saveas('LaMnO₃_pow_sw.sqw','mcstas')



Now is the time to ... mix



Check data files

- Edit the LaMnO3.laz file. Something wrong ?
 - Molar mass: 242 g/mol for LaMnO3
 - $\sigma_{\text{coh}} = 23.0$ barns / LaMnO3
 - $\sigma_{\text{abs}} = 22.3$
 - $\sigma_{\text{inc}} = 1.53$
 - Density = 6.56 g/cm³
- Check other files



Use *LaMnO₃* data files

- We start from the template TOF.
- We remove the container, then add
 - Isotropic_Sqw(Sqw_coh="coh", Sqw_inc="inc", norm=1)
 - Isotropic_Sqw(Sqw_coh="sw", norm=1)
 - PowderN(reflections="laz")
- **Sqw**: **T=-1** (no detailed balance), and **norm=1** (data files as is).
- Use Monitor_Sqw to get (q,w) space



S classical – S quantum

- $S^*(q,w) = S(q,w) e^{-\hbar w/2kT}$
- If $T=10\text{ K} \sim 1\text{ meV}$, and $w \sim 30\text{ meV}$
 - $e^{-\hbar w/2kT} \sim e^{30} = 10^{13} \dots$
- This is solved in iFit, but not in Isotropic_Sqw
 - $T=-1$



RUN!

- templateTOF_LaMnO3.instr
- ncount=1e8. Takes 10 minutes.

