



McStas: Liquids and Powders Dynamics

E. Farhi, *ILL*



Isotropic density samples in McStas

Outline

- Scattering law, a bit of theory
- The *Isotropic_Sqw* component and input data files
- How to get $S(q,w)$ data sets:
 - From MD
 - From lattice dynamics
 - From experiments
- Exercise: a “liquid” TOF spectrometer
- Exercise: separate contributions

Disclaimer: in case of errors and uncertainties, please correct me...






Neutron-matter interaction

Selections rules

$$K_f = K_i + Q$$

$$E_f = E_i + \omega$$



Bragg's law (diffraction on structure – atoms separated by distance d)

$$n\lambda = 2\pi / K_i = 2d \sin \theta$$

Scattering law (intensity per solid angle and energy, dynamics)

Holy Book (Squires)

$$\frac{d^2 \sigma}{d\Omega dE_f} = \frac{K_f}{K_i} \left[\frac{\sigma}{4\pi} S(Q, \omega) \right]$$

Dynamical structure factor $S(Q, \omega)$ is characteristic of each material
Reflects ordering of matter (atom/molecule positions – movements - domains)



Scattering law


The double differential scattering cross section gives the probability for a neutron to scatter for a given solid angle and energy

The total intensity is the full integral over all scattering possibilities.

Effective total scattering
cross section

$$\hat{\sigma} = \iint \frac{d^2\sigma}{d\Omega dE_f} d\Omega dE_f$$

V.F. Sears. Adv. Phys.,
24, 1, 1975.



Computing the total scattering probability



The total scattering cross section is given in (Ω, E_f) space, but S is given in (q, ω)
A variable change must be done for the integration (Jacobian).

We like to play games
in (q, ω) space

$$\frac{d\Omega}{d\theta} = -2\pi \sin\theta$$

$$\frac{dq}{d\theta} = -\frac{k_i k_f \sin\theta}{q}$$

Effective cross section
in (q, ω) space

$$\hat{\sigma} = \sigma \iint \frac{S(q, \omega) q}{2k_i^2} dq d\omega$$



Probability to transmit

$$p = e^{-\rho \hat{\sigma} x}$$

Scattering distribution


$$S(q, \omega)$$

with importance sampling to
scatter preferably where S is large





How to get $S(q,w)$ data sets



MD step is done prior to the virtual experiment (NAMD, VASP, GROMACS, ...).
Computationally intensive (e.g. use clusters). Then use FFT(r,t)

McStas provides a few sample $S(Q,w)$: Rb, Ge, H₂O, D₂O, D₂, ...

Isotropic_Sqw: Handles elastic and inelastic for both coherent and incoherent channels



Isotropic_Sqw syntax

- Isotropic_Sqw(
 Sqw_coh=FILE_COH,
 Sqw_inc=FILE_INC,
 radius=R,
 height=H)
• More component parameters can specify geometry, physical properties, ...
• The data files specify the $S(Q, \omega)$ or $S(Q)$ values as a matrix with Q, ω extent. Additional fields can be included as meta data (# lines).





Isotropic_Sqw data format

- Open McStas, select Help / Comp Ref, then Data files. Select a .sqw file, e.g. He4_liq_coh.sqw
- Look at format: Header and meta data, axes, matrix
- Used # fields:
 - density, weight, sigma, Temperature
- Meta data can be given as component parameters.





Isotropic_Sqw data format

```
# Sqw data file for Isotropic_Sqw
# liquid He4: coherent part, no incoherent, atomic number 2
# Elementary Excitation Data by R.J. Donnelly et al., J. Low Temp. Phys., 44 (1981) 471
# WARNING: line width is constant, intensity is not right
#
# Physical parameters:
# V_rho      0.072   atom density per Angs^3
# weight     4.002   in [g/mol]
# density     0.4784  in [g/cm^3]
# sigma_abs  0.00747 absorption scattering cross section in [barn]
# sigma_coh   1.34    coherent scattering cross section in [barn]
# sigma_inc   0       incoherent scattering cross section in [barn]
# Temperature 2       in [K]
# classical   0       experimental, contains Bose factor
#
# q axis values
# vector of m values in Angstroem-1
0.001000 0.011000 0.02 ...
# w axis values
# vector of n values in meV
0.001391 0.011391 0.021391 0.0313 ...
# sqw values (one line per q axis value)
# matrix of S(q,w) values (m rows x n values), one line per q value
9.721422 10.599145 11.344954 ...
```

He4_liq_coh.sqw



Such files can be written from an iFit iData_Sqw2D object:

- `sqw2d = iData_Sqw2D(...);`
- `saveas(sqw2d, 'filename.sqw', 'mcstas')`



Sqw file from experiments

- Reduce the experimental data, correct it for e.g. empty cell, parallax, detector efficiency, ...
- Correct for absorption, incoherent scattering, multiple-scattering, ...
- Integrate over $|Q|$ as *Isotropic_Sqw* is isotropic.
- Write a text file with $[Q, \omega]$ and $S(Q, \omega)$ vectors/matrices. Add meta-data.





Sqw file from MD

- Set the system box, (pseudo) potential/FF
- Equilibrate, and couple to thermostat (NVT)
- Run calculator in MD mode (NVE)
- Convert trajectory (\mathbf{R}, t) into $S(|Q|, \omega)$ [MDANSE]
 - Analysis/Scattering/DCSF and DISF \rightarrow .nc file
- Export into .sqw text files.
 - `sqw2d = iData_Sqw2D('DCSF.nc');`
 - `saveas(sqw2d, 'DCSF.sqw', 'mcstas');`





Sqw file from Lattice Dynamics

- Create a crystal, calculate displacements or so.
- Assemble the 'phonon' representation
- Evaluate $S(Q,w)$ on powder average
- Export to McStas Sqw






Sqw: Exercise 1: a liquid TOF

Aim: A simple spectrometer (and diffractometer)


- Create a new instrument from 'template (test)'.
- Call it **Liquid_simple** and define input parameters (**lambda**=2.36, **string coh**="Rb_liq_coh.sqw", **string inc**="Rb_liq_inc.sqw")
- Insert a **Source_simple** $\phi 1\text{cm}$ sending $\lambda=\text{lambda}$ with $d\lambda/\lambda=1\%$. Focus onto a $1\text{x}1\text{cm}^2$ area.
- Insert an **Isotropic_Sqw** 3m away, using $\sigma_{\text{coh}}=\text{coh}$, $\sigma_{\text{inc}}=\text{inc}$ with $\phi 1\text{cm} \times 5\text{cm}$.





Sqw: Exercise 1: a liquid TOF


- Add a **Monitor_nD** cylindrical detector $\phi 1\text{m} \times 30\text{cm}$, sensitive to (θ, y) for diffraction, centred on the sample, with 100 bins.
- Add the same, but sensitive to **(angle,energy)** with automatic energy limits.
- Save, run in *Trace 3D* to check geometry.





Sqw: Exercise 1: a liquid TOF

- Run in Simulation/PGPLOT mode with $1e8$ neutron events.
- Plot results !
- Comment on the diffraction pattern and the inelastic one.





Sqw: Exercise 2: contributions

- Insert an instrument variable in the DECLARE block, as 'flag_scat'.

```
DECLARE %{\n  int flag_scat=0;\n}%}
```

- After the AT token of the 'sample', insert an EXTEND block that sets flag_scat to the number of SCATTERED events.

```
EXTEND %{\n  flag_scat=SCATTERED; // nb of scattered events\n}%}
```





Sqw: Exercise 2: contributions


- Duplicate the 2 monitors, and make them sensitive to **multiple-scattering** only. Add a WHEN(flag_scatt) between the COMPONENT and the AT keywords:

COMPONENT ...

WHEN(flag_scatt > 1)

AT ...

- Save and Run !
- How much multiple scattering ratio at $\lambda=2.36 \text{ \AA}$
- Make the sample $\phi 5\text{cm}$ and repeat.





Sqw: Exercise 3: from iFit: build

- McStas can be controlled from within iFit.
- Open Matlab/iFit
- Create the Liquid_simple model with:
 - `model = mccode('Liquid_simple')`
- Plot the geometry with:
 - `plot(model)` % has contextual menus
- Then run it with (default 1e6 event and pars)
 - `data = iData(model, [], nan);`
 - `subplot(model)`


iFit



Sqw: Exercise 3: from iFit: eval



- Specify parameters
 - `data = iData(model, 'lambda=2.36; coh=Cu.laz')`
- Do a scan:
 - `data = iData(model, 'lambda=[1.2 2.4 3.6]')`
- Change neutron events #
 - `model.UserData.options.ncount = 1e7;`





Sqw: Exercise 3: from iFit: optim

- Fix all parameters but lambda, Maximize model value:
 - `mlock(model, 'all'); munlock(model, 'lambda')`
 - `xlim(model, 'lambda', [1 3]); % bounds`
 - `fmax(model, 'lambda=2.36', '', nan)`
 - `fmax(model, 'lambda=2.36', 'OutputFcn=fminplot', nan)`



Sqw: Exercise 3: from iFit: crazy



- You can add two McStas models:
 - `model=mccode('instr1')+mccode('instr2')`
- Edit the instrument and re-compile
 - `edit(model)`
- ...

